



TOMITM
ENVIRONMENTAL SOLUTIONS

STERAMIST[®]
POWERED BY BINARY IONIZATION TECHNOLOGY[®]

A **UNIVERSAL** DISINFECTION SOLUTION
FOR **MULTIPLE** DIFFERENT INDUSTRIES

Innovative disinfection & decontamination technology
and services designed for global infection control

CORPORATE INFORMATION

DISCOVER THE **STERAMIST**[®] FAMILY
OF PRODUCTS AND SERVICES TODAY

TOMIMIST.COM
800.525.1698

TOMI™ Environmental Solutions, Inc.

SteraMist® - Innovating for a Safer World®

TOMI™ Environmental Solutions, Inc. (OTCQB:TOMZ) is a global disinfection & decontamination and infection prevention company, providing environmental solutions for indoor surface disinfection through manufacturing, sales and licensing of its premier Binary Ionization Technology® (BIT™) platform, which was invented under a defense grant in association with the Defense Advanced Research Projects Agency (DARPA) of the U.S. Department of Defense. This state of the art technology produces a fog represented by the SteraMist® brand.

TOMI™ Divisions

TOMI™ focuses primarily on four divisions: Life Sciences, Hospital-Healthcare, the TOMI™ Service Network (TSN™), and Food Safety. Our mission is to help our customers create a healthier world through our range of products and services.




HOSPITAL-HEALTHCARE

- Hospitals
- Medical Clinics
- Surgical Centers
- Rehabilitation
- Senior Living Centers
- Emergency Medical Transport



LIFE SCIENCES

- Biotechnology areas (BSL-3 & BSL-4)
- Pharmaceutical Laboratories
- Vivariums
- Research Facilities



TOMI™ SERVICE NETWORK


- Educational Facilities
- Childcare Facilities
- Residential
- Athletic Equipment & Facilities
- Emergency Response Vehicles



FOOD SAFETY

- Food Processing Plants
- Packaging Plants
- Warehouses/Distribution
- Pre/Post-Harvest Areas
- Food Service Facilities


The iHP™ Process (ionized Hydrogen Peroxide)



7.8% hydrogen peroxide BIT™ Solution converts to iHP™ after passing through a cold plasma arc.




iHP™ is carried throughout the mist, moving like a gas throughout the treated area.



iHP™ damages pathogenic organisms through oxidation of proteins, carbohydrates, and lipids.



Cellular disruptions and/or dysfunctions occur and allows for disinfection/decontamination in the targeted area.



- A. Electrodes
- B. Nozzle
- C. Atmospheric Cold Plasma Arc
- D. ionized Hydrogen Peroxide (iHP™)

The Science

BEHIND STERAMIST®

The atmospheric cold plasma arc converts the H_2O_2 molecules into iHP™. As one of the most powerful oxidizing agents in nature, the iHP™ kills the pathogens achieving high efficacy and leaves behind only oxygen and humidity in treated spaces.

Binary Ionization Technology® (BIT™)

SteraMist® BIT™ Efficacy

As an EPA-registered (No. 90150-2) solution and technology combination, our BIT™ solution is subject to multiple efficacy tests and has been tested to comply with the highest GLP standards at independent laboratories. TOMI™ continues updating its BIT™ label to ensure that we meet and exceed EPA and customer standards for quality and results. SteraMist® is also Health Canada Approved (DIN: 02469448).

- Effective broad-spectrum surface and air disinfectant/decontaminant.
- Effective against bacterial spores and gram negative bacteria, including multiple drug resistant organisms (MDROs).
- Does not require adjustment to room temperature or humidity before application.
- Suitable for frequent/daily use.



TOMI™ has validated the efficacy of this process against a wide range of microorganisms, including:

Geobacillus Stearothermophilus	
Clostridioides difficile	MRSA
Pseudomonas	Staphylococcus
Influenza A Virus (h1n1)	Salmonella
Norovirus	Mold and Mildew

Portable Products

FEATURING A PARTICLE SIZE OF **0.05-3 MICRONS**



STERAMIST® SURFACE UNIT

The **SteraMist® Surface Unit** is a fully portable, fast-acting, hand-held, point-and-spray disinfection system. Whether used for specific service needs or for daily use, the single applicator surface unit enables disinfection of all surfaces - including high touch, sensitive equipment, and electronics.

PRODUCT FEATURES

- Application time of only five seconds per square foot, with a seven-minute contact time
- Disinfection of all surfaces, including high touch, sensitive equipment and electronics
- Reaches surfaces that regular disinfectants and manual cleans cannot reach

SteraMist® Surface Unit | Refer to EPA label no. 90150-2 for more information.



STERAMIST® ENVIRONMENT SYSTEM

The **SteraMist® Environment System** is portable and provides complete room disinfection, deodorization, and mold mitigation using multiple treatment applicators (3) per unit.

PRODUCT FEATURES

- Effective whole complete room treatment in just under 45 minutes* for a room (3,663.7 ft.³/104 m³)
- Scalable for small or larger spaces, with no maximum requirements
- Automated/remote-controlled system with downloadable data sets
- Able to be converted into three hand-held Surface Units, extending disinfection reach

SteraMist® Environment System | Refer to EPA label no. 90150-2 for more information.

SteraMist® Services

SteraMist® iHP® Service

In the event of routine maintenance or to controlling a specific threat in an emergency, TOMI™ offers quick and powerful SteraMist® iHP® Service. Ready at any time, the iHP® service team offers room, equipment, and full facility disinfection and decontamination services.

TOMI™ currently provides routine service to some of the most well-respected facilities in the Life Sciences and Hospital-HealthCare industries. Many pharmaceutical companies worldwide rely on SteraMist® to perform their critical, bi-annual facility-wide shutdown.



TSN™ (TOMI™ Service Network)



TSN™ is comprised of professionals who are exclusively licensed and trained to use the SteraMist® Binary Ionization Technology (BIT™) system and solution. SteraMist® provides an EPA-registered tool to assist in a community's disinfection/decontamination, deodorizing, and mold remediation needs. The TOMI™ Service Network (TSN™) **is an innovator in community-based disinfection, decontamination, odor, and mold remediation.**

FRST™ (Forensic Restoration® Service Team)*



When mitigating and recovering from biological threats and crises, **SteraMist® can help.** Certified by the Global Biorisk Advisory Council (GBAC), a division of ISSA, critical jobs such as trauma cleanup and emergency outbreaks require rapid response, and iHP™ technology delivers disinfection & decontamination with thorough efficacy.

*Forensic Restoration is a property of BioSheen Services, LLC.

TOMI™ Training and Support

When you purchase SteraMist®, you're not only making an investment in the **best disinfection and decontamination technology on the market today**, but you're also invited to enter a fully supported partnership with TOMI™ Environmental Solutions, Inc.

As part of this partnership, you unlock a wealth of additional training and support. Based on your facility's specific needs and your purchased product, a certified SteraMist® Technician will conduct an initial site visit. During this visit, training will be conducted with relevant personnel. Afterwards, our technician will evaluate your unique disinfection and decontamination needs. Based upon the evaluation, customized protocols are recommended for both immediate and ongoing support.



TRAINING PROGRAM INCLUDES:

- Hands-on training with staff for application, care, safety, and troubleshooting. Supervisor training included.
- Ongoing phone and digital consultation with a SteraMist® Trainer at training@tomimist.com.
- Basic troubleshooting technique package available upon request during training.
- Site visits upon request, when necessary.

SUPPORT PROGRAM INCLUDES:

- One year warranty with an option to extend.
- 24/7 equipment technical support (maintenance, troubleshooting, calibration) at repair@tomimist.com.
- 24/7 SteraMist® application support (protocols, proper usage) with a SteraMist® Technician at tomiservice@tomimist.com.
- Ongoing protocol consultations with a SteraMist® Technician.



STERAMIST[®]


POWERED BY BINARY IONIZATION TECHNOLOGY[®]



HOSPITAL-HEALTHCARE



LIFE SCIENCES



TOMI[™] SERVICE NETWORK



FOOD SAFETY

TOMI[™] WORLD HEADQUARTERS

9454 Wilshire Boulevard | Penthouse G-6 | Beverly Hills, CA 90212

Phone: 800.525.1698 | Fax: 412.766.7878

TOMI[™] U.S. CORPORATE HEADQUARTERS

8430 Spires Way | Suite N | Frederick, MD 21701

Phone: 800.525.1698 | Fax: 412.766.7878

TOMIMIST.COM | info@tomimist.com

